

Wispa™ 200 and 300 Series Ceiling Mounted Hoist



Self diagnostic - beeps if there is a problem

Freewheel or powered traverse models

Charge indicator LED and audible warning to indicate low battery levels

Ceiling mounted track or freestanding unit options

▶ Soft start/stop for smooth and controlled transfer on all axes

Aesthetically designed case

▲ Variable speed facility on all axes

+ Weight overload limiting device

★ Emergency stop and emergency lowering feature

Simple to recharge

Carry bar with safety hooks and crossover design

Compatible with all Chiltern Invadex patient handling slings and other manufacturers slings with loops subject to risk assessment by a competent professional

● Clear, easy to use handset

X-Y model option allows complete room coverage.



Chiltern Invadex Ltd

Chiltern House
6 Wedgwood Road
Bicester, Oxfordshire
OX26 4UL

Tel: 01869 246470
Fax: 01869 247214
E-mail: sales@chilterninvadex.co.uk
Web: www.chilterninvadex.co.uk

Overseas Enquiries

Tel: +44 (0) 1869 365545
Fax: +44 (0) 1869 365567
E-mail: export@chilterninvadex.co.uk

Chiltern Invadex (NI) Ltd

5 Balloo Way
Balloo Industrial Estate
Bangor, County Down
BT19 7QZ

Tel: 028 9145 0497
Fax: 028 9127 0046

From Eire

Tel: 048 9145 0497
Fax: 048 9127 0046

Chiltern Invadex Ltd

Unit 16 Evans Business Centre
68-74 Queen Elizabeth Avenue
Hillington Industrial Estate
Glasgow
G52 4NQ

Tel: 0141 880 9986
Fax: 0141 880 9990

Versatile

Chiltern Invadex offer an 'XY' system which can provide complete room coverage. For example, a user could be transferred to any bed in a hospital ward. (see figure a).




Figure a

Moving between rooms is also made simple. The Wispa allows easy transfer, as ceiling track can be fitted in most rooms throughout the house (see figure b).




Figure b

Research and design

Two years of professional research, in conjunction with the University of Birmingham / University of Central England and customer focus panels have created this leading edge design.

Service, maintenance and extended warranty programmes

The european standard EN10535 requires that lifting devices shall be serviced every 12 months*. LOLER regulations require an additional inspection every six months.

A nationwide team of trained engineers can quickly install your equipment and service and maintain it to comply with LOLER regulations on a 365 days/24 hour per day basis.

*Subject to conditions - Further information available upon request

Other products available from Chiltern Invadex include:

Showering and Bathing

Shower Cubicles
Level, Easy and Step-in Showers
Shower Accessories
Shower Seats and Shower Chairs
Shower Trolleys and Stretchers
Walk-in and Easy Transfer Baths.

Moving and Handling

Wheelchairs
Mobile Hoists
Slings
Poolside Hoists.



FM 40009 FM 21085

Bicester is BS EN ISO 9001 approved (FM 40009).

Bangor Northern Ireland is BS EN ISO 9002 approved (FM 21085).

Copyright © Chiltern Invadex Ltd

The design of the product shown on this brochure is vested in Chiltern Invadex Ltd and should not be copied or reproduced in anyway whatsoever without the express permission in writing of Chiltern Invadex Ltd.

As part of our commitment to new product development we reserve the right to alter specifications without prior notice.

This is an original design of Chiltern Invadex Ltd in which Design Right subsists.

Published by: Chiltern Invadex Ltd, 6 Wedgwood Road, Bicester, Oxon, OX26 4UL Registered in England No: 1182024 ISSUE 6 Oct 2005

In Northern Ireland: Chiltern Invadex (NI) Ltd, 5 Balloo Way, Balloo Industrial Estate, Bangor, County Down, BT19 7QZ Registered No: NI 19958